



## Anti-Mouse KIRREL3/NEPH2 Monoclonal Antibody, N321C/49R

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19248
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a, kappa
<b>TARGET</b>	KIRREL3, NEPH2, Kin of irregular chiasm-like protein 3, Nephrin-like protein 2, KIAA1867, Kin of IRRE-like protein 3

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Mouse, Rat
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

