



Anti-Human Gal-10/CLC Monoclonal Antibody, 1D11

OVERVIEW

PRODUCT REFERENCE	PTX19254
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Lysolecithin acylhydrolase, Galectin-10, LGALS10A, CLC, Gal-10, LGALS10, Eosinophil lysophospholipase, Charcot-Leyden crystal protein

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

