



Anti-Rat GRIA1/GluR-1 Monoclonal Antibody, 11B8

OVERVIEW

PRODUCT REFERENCE	PTX19255
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Glutamate receptor 1, GluR-1, AMPA-selective glutamate receptor 1, GluR-A, GluR-K1, Glutamate receptor ionotropic, AMPA 1 (GluA1)

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Rat
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

