



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Rat GRIA2/GluR-2 Monoclonal Antibody, 15F1

OVERVIEW

PRODUCT REFERENCE	PTX19256
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Glutamate receptor 2; GluR-2; AMPA-selective glutamate receptor 2; GluR-B; GluR-K2; Glutamate receptor ionotropic, AMPA 2; GluA2; Gria2; Glur2

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Rat
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.