



Anti-Human GluA2 Monoclonal Antibody, L21/32R

OVERVIEW

PRODUCT REFERENCE	PTX19258
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	GluR-2, AMPA-selective glutamate receptor 2, GluR-B, Glutamate receptor 2, GLUR2, GRIA2, GluA2, Glutamate receptor ionotropic, AMPA 2, GluR-K2

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human, Mouse, Rat
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

