



Anti-RSV F/Fusion glycoprotein F0 Monoclonal Antibody, 4D7

OVERVIEW

PRODUCT REFERENCE	PTX19259
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Fusion glycoprotein F0; Fusion glycoprotein F2; F2; p27; Intervening segment; Pep27; Peptide 27; Fusion glycoprotein F1; F1; F; post-fusion

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human respiratory syncytial virus A (strain A2)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

