



Anti-Human VAP-1/AOC3 Monoclonal Antibody, TK8-14

OVERVIEW

PRODUCT REFERENCE	PTX19260
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Membrane primary amine oxidase, AOC3, SSAO, Semicarbazide-sensitive amine oxidase, Vascular adhesion protein 1, VAP-1, HPAO, Copper amine oxidase, VAP1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

