



Anti-Human MAPT/Tau/PHF-tau Monoclonal Antibody, RB86

OVERVIEW

PRODUCT REFERENCE	PTX19263
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Microtubule-associated protein tau, Neurofibrillary tangle protein, Paired helical filament-tau, PHF-tau, MAPT, MAPTL, MTBT1, TAU

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

