



## Anti-FMDV Capsid protein VP1 Monoclonal Antibody, 4C4

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19269
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, kappa
<b>TARGET</b>	Genome polyprotein, Leader protease, Lpro, 3.4.22.46, Capsid protein VP1, P1D, Virion protein 1, Foot-and-mouth disease virus(FMDV)

### PRODUCT INFORMATION

---

<b>HOST</b>	Human
<b>REACTIVITY</b>	Foot-and-mouth disease virus (isolate -/Germany/C1Oberbayen/1960 serotypeC) (FMDV)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

