



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-FMDV Capsid protein VP1 Monoclonal Antibody, 1E12#

OVERVIEW

PRODUCT REFERENCE	PTX19271
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Genome polyprotein, Leader protease, Lpro, 3.4.22.46, Capsid protein VP1, P1D, Virion protein 1, Foot-and-mouth disease virus(FMDV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Foot-and-mouth disease virus (isolate Bovine/Brazil/A24Cruzeiro/1955serotype A) (FMDV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.