



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-AAV2 VP1/Capsid protein VP1 Monoclonal Antibody, SAA0358

OVERVIEW

PRODUCT REFERENCE	PTX19273
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Capsid protein VP1, VP1, Human Adenovirus (HAdV), Adeno-associated virus 2, AAV2, AAV-2

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Adeno-associated virus 2 (isolate Srivastava/1982) (AAV-2)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.