



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human parvovirus B19/HPV-B19 VP2/Viral protein 2 Monoclonal Antibody, SAA0341

OVERVIEW

PRODUCT REFERENCE	PTX19274
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	VP2 protein, VP2, Viral protein 2, Major capsid protein VP2, Human Parvovirus B19 (B19V)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Human parvovirus B19 (HPV B19)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.