



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Coxsackievirus A6/CVA6 VP1/Capsid protein VP1 Monoclonal Antibody, SAA0354

OVERVIEW

PRODUCT REFERENCE	PTX19279
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Genome polyprotein, Capsid protein VP1, P1D, Virion protein 1, Coxsackievirus A6, CVA6

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Coxsackievirus A6
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.