



Anti-Mouse GIPR/GIP-R Monoclonal Antibody, SAA0400

OVERVIEW

PRODUCT REFERENCE	PTX19281
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Gastric inhibitory polypeptide receptor, GIP-R, Glucose-dependent insulinotropic polypeptide receptor, Gipr

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

