



Anti-Human FGF12 Monoclonal Antibody, N94/17R

OVERVIEW

PRODUCT REFERENCE	PTX19285
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Fibroblast growth factor homologous factor 1, Myocyte-activating factor, FHF-1, FGF12B, FHF1, FGF-12, FGF12, Fibroblast growth factor 12

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

