



Anti-Human HRAS Monoclonal Antibody, SAA0505

OVERVIEW

PRODUCT REFERENCE	PTX19292
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	c-H-ras, GTPase HRas, H-Ras-1, p21ras, Transforming protein p21, Ha-Ras, HRAS1, HRAS

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

