



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Anti-Staphylococcus aureus LukAB/LukGH Monoclonal Antibody, SAA0520

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19298
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	Cytotoxin LukAB, Leukocidin LukGH, LukG, LukA, Leukocidin G, Bi-component leukocidin LukGH subunit G, LukH, LukB, Leukocidin H, Bi-component leukocidin LukGH subunit H, Staphylococcus aureus

### PRODUCT INFORMATION

---

<b>HOST</b>	Human
<b>REACTIVITY</b>	Staphylococcus aureus (strain COL), Staphylococcus aureus (strain NCTC 8325 / PS 47), Staphylococcus aureus (strain NCTC 8325 / PS 47), Staphylococcus aureus (strain COL)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**