



Anti-Human oxLDL Monoclonal Antibody, SAA0552

OVERVIEW

PRODUCT REFERENCE	PTX19303
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	oxidised LDL, oxidised form of low-density lipoprotein (LDL), oxidised low-density lipoproteins, oxidized low-density lipoprotein, oxidized low-density lipoproteins, OxLDL

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

