



Anti-Dabigatran Etxilate Monoclonal Antibody, SAA0560

OVERVIEW

PRODUCT REFERENCE	PTX19306
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	dabigatran etxilate, Dabigatran etxilate mesylate, Dabigatran etxilate methanesulfonate

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	General
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

