



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-pan-Aflatoxin Monoclonal Antibody, SAA0544

OVERVIEW

PRODUCT REFERENCE	PTX19312
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	AFM1, AFB1, AFB2, AFG1, AFG2, Aflatoxin M1, 4-Hydroxyaflatoxin B1, Aflatoxin B1, Aflatoxin G1, Aflatoxin G2, Aflatoxin B2, Aflatoxin

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	AFM1, AFB1, AFB2, AFG1, AFG2
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.