



## Anti-Human NECTIN4 Monoclonal Antibody, P1E7

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19319
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1, lambda
<b>TARGET</b>	LNIR, NECTIN4, PRR4, Poliovirus receptor-related protein 4, Nectin cell adhesion molecule 4, Ig superfamily receptor LNIR, Nectin-4, PVRL4

### PRODUCT INFORMATION

---

<b>HOST</b>	Human
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

