



Anti-Mouse NECTIN4 Monoclonal Antibody, R4P1

OVERVIEW

PRODUCT REFERENCE	PTX19320
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, lambda
TARGET	LNIR, NECTIN4, PRR4, Poliovirus receptor-related protein 4, Nectin cell adhesion molecule 4, Ig superfamily receptor LNIR, Nectin-4, PVRL4

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

