



Anti-Human SULT1E1 Monoclonal Antibody, SAA0758

OVERVIEW

PRODUCT REFERENCE	PTX19321
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG
TARGET	ST1E1, STE, EST-1, SULT1E1, Sulfotransferase 1E1, Sulfotransferase, estrogen-preferring, Estrogen sulfotransferase

PRODUCT INFORMATION

HOST	Rabbit
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

