



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human TRPM4/LTRPC4 Monoclonal Antibody, SAA0760

OVERVIEW

PRODUCT REFERENCE	PTX19323
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-hFc
TARGET	LTRPC4, Transient receptor potential cation channel subfamily M member 4, TRPM4, Melastatin-4, LTrpC-4, hTRPM4, Calcium-activated non-selective cation channel 1, LTrpC4, Long transient receptor potential channel 4

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.