



## Anti-Human SERPINE1/PAI-1 Monoclonal Antibody, 8H9D4

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19324
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a
<b>TARGET</b>	PLANH1, Endothelial plasminogen activator inhibitor, SERPINE1, Plasminogen activator inhibitor 1, PAI1, PAI, PAI-1, Serpin E1

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human, Pig
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

