



## Anti-Clostridium tetani tetX/Tetanus toxin Monoclonal Antibody, SAA0801

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19325
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a
<b>TARGET</b>	Tetanus toxin, 3.4.24.68, Tentyoxylisin, Tetanus toxin light chain, Tetanus toxin chain L, Tetanus toxin heavy chain, Tetanus toxin chain H, tetX, CTC_p60

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Clostridium tetani (strain Massachusetts / E88)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

