



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human ADGRE1/EMR1 Monoclonal Antibody, SAA0811

OVERVIEW

PRODUCT REFERENCE	PTX19327
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	ADGRE1, EGF-like module-containing mucin-like hormone receptor-like 1, EGF-like module receptor 1, Adhesion G protein-coupled receptor E1, EMR1, TM7LN3, EMR1 hormone receptor

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.