



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Vibrio cholerae ctxB/Cholera Toxin Subunit B Monoclonal Antibody, SAA0849

OVERVIEW

PRODUCT REFERENCE	PTX19333
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Cholera enterotoxin B-subunit, Cholera enterotoxin subunit B, Cholera toxin B protein (CTB), Cholera toxin B subunit, Cholera toxin beta subunit, Cholera toxin subunit B, Cholerae toxin B subunit, CtxB

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Vibrio cholerae
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.