



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Clostridium tetani tetX/Tetanus toxin Monoclonal Antibody, SAA0858

OVERVIEW

PRODUCT REFERENCE	PTX19339
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Tetanus toxin, 3.4.24.68, Tentyoxylisin, Tetanus toxin light chain, Tetanus toxin chain L, Tetanus toxin heavy chain, Tetanus toxin chain H, tetX, CTC_p60

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Clostridium tetani (strain Massachusetts / E88)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.