



Anti-Clostridium tetani tetX/Tetanus toxin Monoclonal Antibody, SAA0860

OVERVIEW

PRODUCT REFERENCE	PTX19341
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Tetanus toxin, 3.4.24.68, Tentyoxylisin, Tetanus toxin light chain, Tetanus toxin chain L, Tetanus toxin heavy chain, Tetanus toxin chain H, tetX, CTC_p60

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Clostridium tetani (strain Massachusetts / E88)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

