



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-CPV-2 Capsid protein VP1 Monoclonal Antibody, SAA0843

OVERVIEW

PRODUCT REFERENCE	PTX19349
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	Capsid protein VP1, Coat protein VP1, Canine parvovirus type 2 (CPV-2)

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Canine parvovirus type 2 (isolate Dog/United States/CPV-d/1988) (CPV-2), Feline panleukopenia virus (FPV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.