



## Anti-Human S100A9 Monoclonal Antibody, 6B4

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19350
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a, kappa
<b>TARGET</b>	Protein S100-A9, CAGB, MRP-14, S100 calcium-binding protein A9, Leukocyte L1 complex heavy chain, Migration inhibitory factor-related protein 14, p14, Calprotectin L1H subunit, S100A9, Calgranulin-B, CFAG, MRP14

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

