



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Influenza A virus M/Matrix protein 2 Monoclonal Antibody, 14C2

OVERVIEW

PRODUCT REFERENCE	PTX19355
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
TARGET	Matrix protein 2, Proton channel protein M2, M, Influenza Virus (Flu)

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Influenza A virus (strain A/Puerto Rico/8/1934 H1N1)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.