



Anti-Bovine Serum Albumin/BSA Monoclonal Antibody, B1-4

OVERVIEW

PRODUCT REFERENCE	PTX19358
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Albumin, BSA, Bos d 6, ALB

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Bovine
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

