



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Human CD124/IL4R/IL-4R α Monoclonal Antibody, SAA1420

OVERVIEW

PRODUCT REFERENCE	PTX19359
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	IL-4RA, IL-4-binding protein, IL-4 receptor subunit alpha, Soluble IL-4R-alpha, IL4-BP, IL4R, IL-4R subunit alpha, IL-4R-alpha, CD124, Interleukin-4 receptor subunit alpha, IL4RA, Soluble IL-4 receptor subunit alpha, sIL4Ralpha/prot

PRODUCT INFORMATION

HOST	
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.