



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Pan HBV Pres1/Pre-S1 Region Monoclonal Antibody, KR127

OVERVIEW

PRODUCT REFERENCE	PTX19361
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Large envelope protein S, L glycoprotein, L-HBsAg, LHB, Large surface protein, Major surface protein, Hepatitis B Virus (HBV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Hepatitis B virus (HBV) (Q2EID8), Hepatitis B virus genotype A2 subtype adw (isolate Japan/Nishioka/1983)(HBV-A) (P03142), Hepatitis B virus genotype D subtype adw (isolate United Kingdom/adyw/1979)(HBV-D) (P03139), Hepatitis B virus (HBV) (Q96836), Hepa
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.