



Anti-HBV-D HBeAg/C/p25 Monoclonal Antibody, E6#

OVERVIEW

PRODUCT REFERENCE	PTX19363
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	External core antigen, HBeAg, Precore protein, p25, C, Hepatitis B Virus (HBV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Hepatitis B virus genotype D subtype ayw (isolate France/Tiollais/1979) (HBV-D)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

