



Anti-HBV-D HBcAg/HBeAg Monoclonal Antibody, E13#

OVERVIEW

PRODUCT REFERENCE	PTX19364
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Capsid protein, Core antigen, Core protein, HBcAg, p21.5, C, Hepatitis B Virus (HBV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Hepatitis B virus genotype D subtype adw (isolate United Kingdom/adyw/1979)(HBV-D) (P03147), Hepatitis B virus genotype D subtype ayw (isolate France/Tiollais/1979)(HBV-D)(P0C573)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

