



Anti-HIV-1 Capsid protein p24/CA protein Monoclonal Antibody, A3

OVERVIEW

PRODUCT REFERENCE	PTX19366
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Pr160Gag-Pol, Gag-Pol polyprotein, gag-pol, Capsid protein p24, CA, Human Immunodeficiency Virus (HIV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Human immunodeficiency virus type 1 group M subtype B (isolate BH10) (HIV-1)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

