



## Anti-HIV-1 rev/ART/TRS protein Monoclonal Antibody, SAA1422

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19370
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	Protein Rev, ART/TRS, Anti-repression transactivator, Regulator of expression of viral proteins, rev, Human Immunodeficiency Virus (HIV)

### PRODUCT INFORMATION

---

<b>HOST</b>	
<b>REACTIVITY</b>	Human immunodeficiency virus type 1 group M subtype B (isolate HXB2) (HIV-1)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

