



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-HIV-1 nef/F-protein protein Monoclonal Antibody, SAA1423

OVERVIEW

PRODUCT REFERENCE	PTX19371
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	Protein Nef, 3'ORF, Negative factor, F-protein, C-terminal core protein, nef, Human Immunodeficiency Virus (HIV)

PRODUCT INFORMATION

HOST	
REACTIVITY	Human immunodeficiency virus type 1 group M subtype B (isolate ARV2/SF2) (HIV-1)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.