



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Influenza A virus NP/Nucleoprotein protein Monoclonal Antibody, SAA1425

OVERVIEW

PRODUCT REFERENCE	PTX19372
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	Nucleoprotein, Nucleocapsid protein, Protein N, NP, Influenza Virus (Flu)

PRODUCT INFORMATION

HOST	
REACTIVITY	Influenza A virus (strain A/Wilson-Smith/1933 H1N1) (Influenza A virus(strain A/WS/1933 H1N1)) (Q1K9H2), Influenza A virus (A/WSN/1933(H1N1)) (Q9DLK6)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.