



## Anti-Influenza A virus RNA polymerase complex Monoclonal Antibody, SAA1427

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19373
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	Polymerase acidic protein, 3.1.-., RNA-directed RNA polymerase subunit P2, PA, Influenza Virus (Flu)

### PRODUCT INFORMATION

---

<b>HOST</b>	
<b>REACTIVITY</b>	Influenza A virus (strain A/Northern Territory/60/1968 H3N2) (Influenza A virus (strain NT60)) (Influenza A virus (strain A/NT/60/1968 H3N2)) (P03434 & P03432), Influenza A virus (A/Northern Territory/60/1968(H3N2)) (Q1AR05 & Q910D6)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

