



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Candida albicans ENO1/Enolase 1 Monoclonal Antibody, SAA1432

OVERVIEW

PRODUCT REFERENCE	PTX19374
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Enolase 1, 4.2.1.11, 2-phospho-D-glycerate hydro-lyase, 2-phosphoglycerate dehydratase, ENO1, CAALFM_C108500CA, CaO19.395, CaO19.8025

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast) (P30575), Candida albicans (Yeast) (V5UZY2), Candida tropicalis (strain ATCC MYA-3404 / T1) (Yeast) (C5MAS1)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.