



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Mouse CD223/LAG3 Monoclonal Antibody, SAA0767

OVERVIEW

PRODUCT REFERENCE	PTX19375
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	FDC, LAG3, Lymphocyte activation gene 3 protein, CD223, sLAG-3, LAG-3, FMD30410

PRODUCT INFORMATION

HOST	
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.