



Anti-Ganglioside GD2 Monoclonal Antibody, ME36.1

OVERVIEW

PRODUCT REFERENCE	PTX19377
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	Ganglioside GD2, CAS: 65988-71-8

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	General
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

