



Anti-Human CACNA1C Monoclonal Antibody, N263/31R

OVERVIEW

PRODUCT REFERENCE	PTX19380
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	Cacna1c, Cach2, Cacn2, Cacn1a1, Cchl1a1, Voltage-dependent L-type calcium channel subunit alpha-1C, Calcium channel, L type, alpha-1 polypeptide, isoform 1, cardiac muscle, Rat brain class C, RBC, Voltage-gated calcium channel subunit alpha Cav1.2

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

