



## Anti-Rabbit CACNA1C Monoclonal Antibody, L57/23R

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19382
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a
<b>TARGET</b>	Target name:CACNA1C, CACH2, CACN2, CACNL1A1, CCHL1A1, Voltage-dependent L-type calcium channel subunit alpha-1C, Calcium channel L type alpha-1 polypeptide, Smooth muscle calcium channel blocker receptor, CACB-receptor, Voltage-gated calcium channel subun

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Rabbit
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

