



Anti-Human TRIM33 Monoclonal Antibody, SAA1438

OVERVIEW

PRODUCT REFERENCE	PTX19383
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a
TARGET	TRIM33, KIAA1113, RFG7, TIF1G, E3 ubiquitin-protein ligase TRIM33, Ectodermin homolog, RET-fused gene 7 protein, Protein Rfg7, RING-type E3 ubiquitin transferase TRIM33, Transcription intermediary factor 1-gamma, TIF1-gamma, Tripartite motif-containing pr

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

