



Anti-Mouse MAPT/Tau/PHF-tau Monoclonal Antibody, TAU-5

OVERVIEW

PRODUCT REFERENCE	PTX19384
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	MAPTL, Paired helical filament-tau, MTBT1, Microtubule-associated protein tau, TAU, Neurofibrillary tangle protein, MAPT, PHF-tau

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Mouse
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

